

## Product Datasheet

### Human PAPP A ELISA Kit (orb406985)

<b>Description</b>	Human pregnancy associated plasma protein A,PAPP-A ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.7 ng/mL-50 ng/mL
<b>Target</b>	PAPP A
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, urine
<b>Uniprot ID</b>	<a href="#">Q13219</a>
<b>Sensitivity</b>	0.4 ng/mL
<b>Expiration Date</b>	Please enquire.