

## Product Datasheet

### Human TFF2 ELISA Kit (orb407010)

<b>Description</b>	Human Trefoil factor 2(TFF2) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 ng/mL-2000 ng/mL
<b>Target</b>	TFF2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernates
<b>Uniprot ID</b>	<b>Q03403</b>
<b>Sensitivity</b>	7.81 ng/mL
<b>Expiration Date</b>	Please enquire.