

Product Datasheet

Human DHH ELISA Kit (orb407017)

| | |
|--------------------------|---|
| Description | Human Desert hedgehog (DHH) ELISA kit |
| Reactivity | Human |
| Range | 15.6 pg/mL-1000 pg/mL |
| Target | DHH |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, cell culture supernates, tissue homogenates |
| Uniprot ID | O43323 |
| Sensitivity | 3.9 pg/mL |
| Expiration Date | Please enquire. |