

Product Datasheet

Human CXCL1 ELISA Kit (orb407074)

Description	Human growth-regulated oncogene α /melanoma growth stimulating activity,GRO α /MGSA ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Target	CXCL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09341
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.