

## Product Datasheet

### Human GAS6 ELISA Kit (orb407076)

**Description** Human Growth arrest-specific protein 6(GAS6) ELISA kit

**Reactivity** Human

**Range** 0.78 ng/mL-50 ng/mL

**Target** GAS6

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [Q14393](#)

**Sensitivity** 0.195 ng/mL

**Expiration Date** Please enquire.