

Product Datasheet

Human CHGA ELISA Kit (orb407082)

Description	Human Chromogranin,CgA ELISA Kit
Reactivity	Human
Range	62.5 ng/mL-4000 ng/mL
Target	CHGA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, cerebrospinal fluid (CSF)
Uniprot ID	P10645
Sensitivity	15.6 ng/mL
Expiration Date	Please enquire.