

## Product Datasheet

### Human PPIB ELISA Kit (orb407083)

**Description** Human cyclophilin B,CyPB ELISA Kit

**Reactivity** Human

**Range** 31.25 pg/mL-2000 pg/mL

**Target** PPIB

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell lysates, tissue homogenates

**Uniprot ID** [P23284](#)

**Sensitivity** 7.8 pg/mL

**Expiration Date** Please enquire.