

Product Datasheet

Human Aquaporin 1 ELISA Kit (orb407093)

| | |
|--------------------------|--|
| Description | Human Aquaporin 1 ELISA Kit |
| Reactivity | Human |
| Range | 15.6 pg/mL-1000 pg/mL |
| Target | AQP1 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, cell culture supernates, cell lysates |
| Uniprot ID | P29972 |
| Sensitivity | 3.9 pg/mL |
| Expiration Date | Please enquire. |