

Product Datasheet

Human Aquaporin 4 ELISA Kit (orb407096)

| | |
|--------------------------|-----------------------------------|
| Description | Human Aquaporin 4 ELISA Kit |
| Reactivity | Human |
| Range | 0.156 ng/mL-10 ng/mL |
| Target | AQP4 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P55087 |
| Sensitivity | 0.039 ng/mL |
| Expiration Date | Please enquire. |