

## Product Datasheet

### Human RLN1 ELISA Kit (orb407118)

<b>Description</b>	Human Relaxin,RLN ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	6.25 pg/mL-400 pg/mL
<b>Target</b>	RLN1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P04808</b>
<b>Sensitivity</b>	1.56 pg/mL
<b>Expiration Date</b>	Please enquire.