

## Product Datasheet

### Human TREM1 ELISA Kit (orb407127)

<b>Description</b>	Human Triggering Receptor Expresses on Myeloid Cells-1,TREM-1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	TREM1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9NP99</a>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.