

## Product Datasheet

### Human PGF ELISA Kit (orb407138)

<b>Description</b>	Human placenta growth factor,PLGF ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	1.67 ng/mL-40 ng/mL
<b>Target</b>	PGF
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P49763</b>
<b>Sensitivity</b>	0.7 ng/mL
<b>Expiration Date</b>	Please enquire.