

## Product Datasheet

### Human GOT1 ELISA Kit (orb407166)

<b>Description</b>	Human Aspartate aminotransferase,AST ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312 mU/ml - 20 mU/ml
<b>Target</b>	GOT1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P17174</b>
<b>Sensitivity</b>	0.078 mU/ml
<b>Expiration Date</b>	Please enquire.