

## Product Datasheet

### Human AZU1 ELISA Kit (orb407168)

<b>Description</b>	Human Azurocidin/Heparin Binding Protein(AZU/HBP) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8 pg/mL-500 pg/mL
<b>Target</b>	AZU1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, urine, tissue homogenates, cell culture supernates
<b>Uniprot ID</b>	<b>P20160</b>
<b>Sensitivity</b>	1.95 pg/mL
<b>Expiration Date</b>	Please enquire.