

## Product Datasheet

### Human Ferritin ELISA Kit (orb407170)

<b>Description</b>	Human Ferritin, Heavy Polypeptide (FTH) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.45 ng/mL-30 ng/mL
<b>Target</b>	FTH1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P02794</b>
<b>Sensitivity</b>	0.11 ng/mL
<b>Expiration Date</b>	Please enquire.