

Product Datasheet

Human Tubulin beta 3 ELISA Kit (orb407204)

Description	Human Tubulin Beta-3 Chain (TUBB3)ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	TUBB3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13509
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.