

Product Datasheet

Human MAP1LC3B ELISA Kit (orb407206)

Description	Human Microtubule-associated proteins 1A/1B light chain 3B(MAP1LC3B) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Target	MAP1LC3B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9GZQ8
Sensitivity	4.67 pg/mL
Expiration Date	Please enquire.