

Product Datasheet

Human ProGRP ELISA Kit (orb407220)

| | |
|--------------------------|---|
| Description | Human pro-gastrin-releasing peptide, ProGRP ELISA Kit |
| Reactivity | Human |
| Range | 18.75 pg/mL-1200 pg/mL |
| Target | ProGRP |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Sensitivity | 4.68 pg/mL |
| Expiration Date | Please enquire. |