

## Product Datasheet

### Human cellular fibronectin ELISA Kit (orb407261)

<b>Description</b>	Human cellular fibronectin,cFn ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	12.5 ng/mL-800 ng/mL
<b>Target</b>	cFn
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<a href="#">P02751</a>
<b>Sensitivity</b>	3.12 ng/mL
<b>Expiration Date</b>	Please enquire.