

Product Datasheet

Human AFP-L3 ELISA Kit (orb407285)

Description	Human alpha-fetoprotein Lens culinaris agglutinin 3,AFP-L3 ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	AFP-L3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Sensitivity	0.3 ng/mL
Expiration Date	Please enquire.