

## Product Datasheet

### Human CTF1 ELISA Kit (orb407294)

**Description** Human cardiostrophin 1, CT-1 ELISA Kit

**Reactivity** Human

**Range** 7.8 pg/mL-500 pg/mL

**Target** CTF1

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell culture supernates,  
tissue homogenates

**Uniprot ID** [Q16619](#)

**Sensitivity** 1.95 pg/mL

**Expiration Date** Please enquire.