

Product Datasheet

Human FLT1 ELISA Kit (orb407331)

| | |
|--------------------------|---|
| Description | Human vascular endothelial cell growth factor receptor 1 (VEGFR-1/Flt1) ELISA kit |
| Reactivity | Human |
| Range | 0.156 ng/mL-10 ng/mL |
| Target | VEGFR-1/Flt1 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P17948 |
| Sensitivity | 0.039 ng/mL |
| Expiration Date | Please enquire. |