

## Product Datasheet

### Human VEGFC ELISA Kit (orb407333)

**Description** Human Vascular Endothelial cell Growth Factor C, VEGF-C ELISA KIT

**Reactivity** Human

**Range** 15.6 pg/mL-1000 pg/mL

**Target** VEGFC

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P49767](#)

**Sensitivity** 3.9 pg/mL

**Expiration Date** Please enquire.