

Product Datasheet

Human ANGPTL2 ELISA Kit (orb407341)

Description	Human angiotensin-like protein 2(ANGPTL2) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	ANGPTL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9UKU9
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.