

Product Datasheet

Human ANGPTL3 ELISA Kit (orb407342)

Description	Human angiotensin-like protein 3, ANGPTL3 ELISA Kit
Reactivity	Human
Range	78 pg/mL-5000 pg/mL
Target	ANGPTL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Y5C1
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.