

## Product Datasheet

### Human PRDX4 ELISA Kit (orb407395)

**Description** Human Peroxiredoxin 4(PRDX4)  
ELISA Kit

**Reactivity** Human

**Range** 0.625 ng/mL-40 ng/mL

**Target** PRDX4

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma

**Uniprot ID** [Q13162](#)

**Sensitivity** 0.156 ng/mL

**Expiration Date** Please enquire.