

## Product Datasheet

### Human IGF-1 ELISA Kit (orb407409)

<b>Description</b>	Human insulin-like growth factors 1,IGF-1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8 ng/mL-500 ng/mL
<b>Target</b>	IGF-1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P05019</b>
<b>Sensitivity</b>	1.95 ng/mL
<b>Expiration Date</b>	Please enquire.