

Product Datasheet

Human IGF2BP2 ELISA Kit (orb407411)

Description	Human Insulin-like growth factor 2 mRNA-binding protein 2(IGF2BP2) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	IGF2BP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y6M1
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.