

## Product Datasheet

### Human IGFBP1 ELISA Kit (orb407415)

<b>Description</b>	Human insulin-like growth factors binding protein 1,IGFBP-1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56 ng/mL-100 ng/mL
<b>Target</b>	IGFBP1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P08833</b>
<b>Sensitivity</b>	0.39 ng/mL
<b>Expiration Date</b>	Please enquire.