

Product Datasheet

Human IGFBP4 ELISA Kit (orb407418)

Description	Human insulin-like growth factors binding protein 4,IGFBP-4 ELISA Kit
Reactivity	Human
Range	62.5 ng/mL-4000 ng/mL
Target	IGFBP4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P22692
Sensitivity	15.6 ng/mL
Expiration Date	Please enquire.