

## Product Datasheet

### Human IGFBP7 ELISA Kit (orb407420)

**Description** Human Insulin-like growth factor-binding protein 7 (IGFBP7/MAC25/PSF) ELISA kit

**Reactivity** Human

**Range** 78 pg/ml - 5000 pg/ml

**Target** IGFBP7

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, homogenates

**Uniprot ID** [Q16270](#)

**Sensitivity** 19.5 pg/ml

**Expiration Date** Please enquire.