

Product Datasheet

Human MT-RNR2 ELISA Kit (orb407470)

Description	Human Putative humanin peptide(MT-RNR2) ELISA kit
Reactivity	Human
Range	28 pg/mL-1800 pg/mL
Target	MT-RNR2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8IVG9
Sensitivity	7 pg/mL
Expiration Date	Please enquire.