

## Product Datasheet

### Human CCL28 ELISA Kit (orb407511)

<b>Description</b>	Human mucosae associated epithelia chemokine,MEC ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78 pg/ml - 5000 pg/ml
<b>Target</b>	CCL28
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>Q9NRJ3</b>
<b>Sensitivity</b>	19.5 pg/ml
<b>Expiration Date</b>	Please enquire.