

## Product Datasheet

### Human TNFRSF10B ELISA Kit (orb407554)

<b>Description</b>	Human Tumor Necrosis Factor-Related Apoptosis-Inducing Ligand Receptor 2 (TRAIL-R2/DR5) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	TNFRSF10B
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Uniprot ID</b>	<a href="#">O14763</a>
<b>Sensitivity</b>	7.8 pg/mL
<b>Expiration Date</b>	Please enquire.