

## Product Datasheet

### Human TNFRSF10D ELISA Kit (orb407556)

|                          |                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human tumor necrosis factor-related apoptosis-inducing ligand receptor 4 (TRAIL-R4) ELISA Kit |
| <b>Reactivity</b>        | Human                                                                                         |
| <b>Range</b>             | 0.156 ng/ml - 10 ng/ml                                                                        |
| <b>Target</b>            | TNFRSF10D                                                                                     |
| <b>Note</b>              | For research use only                                                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                                             |
| <b>Uniprot ID</b>        | <b>Q9UBN6</b>                                                                                 |
| <b>Sensitivity</b>       | 0.039 ng/ml                                                                                   |
| <b>Expiration Date</b>   | Please enquire.                                                                               |