

## Product Datasheet

### Human NTN1 ELISA Kit (orb407562)

<b>Description</b>	Human Netrin-1,Ntn1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	NTN1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Uniprot ID</b>	<b>O95631</b>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.