

Product Datasheet

Mouse ANGPTL2 ELISA Kit (orb407758)

Description	Mouse Angiopoietin-related protein 2(ANGPTL2) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	ANGPTL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9R045
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.