

## Product Datasheet

### Mouse BTLA ELISA Kit (orb407787)

<b>Description</b>	Mouse B- and T-lymphocyte attenuator(BTLA) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.625 ng/mL-40 ng/mL
<b>Target</b>	BTLA
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q7TSA3</a>
<b>Sensitivity</b>	0.156 ng/mL
<b>Expiration Date</b>	Please enquire.