

Product Datasheet

Mouse Interferon gamma ELISA Kit (orb407968)

Description	Mouse C-X-C motif chemokine 9(CXCL9) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	CXCL9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P18340
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.