

## Product Datasheet

### Mouse EFEMP1 ELISA Kit (orb408031)

<b>Description</b>	Mouse EGF-containing fibulin-like extracellular matrix protein 1(EFEMP1) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	EFEMP1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<a href="#">Q8BPB5</a>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.