

Product Datasheet

Mouse EFNA1 ELISA Kit (orb408055)

| | |
|--------------------------|-----------------------------------|
| Description | Mouse Ephrin-A1(EFNA1) ELISA kit |
| Reactivity | Mouse |
| Range | 93.75 pg/mL-6000 pg/mL |
| Target | EFNA1 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P52793 |
| Sensitivity | 23.43 pg/mL |
| Expiration Date | Please enquire. |