

Product Datasheet

Mouse EGFL7 ELISA Kit (orb408057)

| | |
|--------------------------|---|
| Description | Mouse Epidermal growth factor-like protein 7(EGFL7) ELISA kit |
| Reactivity | Mouse |
| Range | 7.8 ng/mL-500 ng/mL |
| Target | EGFL7 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma,tissue homogenates |
| Uniprot ID | Q9QXT5 |
| Sensitivity | 1.95 ng/mL |
| Expiration Date | Please enquire. |