

## Product Datasheet

### Mouse IL22RA2 ELISA Kit (orb408261)

<b>Description</b>	Mouse Interleukin-22 receptor subunit alpha-2(IL22RA2) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	46.88 pg/mL-3000 pg/mL
<b>Target</b>	IL22RA2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q80XF5</a>
<b>Sensitivity</b>	11.7 pg/mL
<b>Expiration Date</b>	Please enquire.