


Product Datasheet

Cystatin C antibody (orb40828)

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Cystatin C antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, IP, WB |
| Immunogen | Recombinant Human Troponin I, fast skeletal muscle protein (1-100AA) |
| Target | TNNI2 |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P48788 |
| Dilution Range | IHC-P:1:20-200, IP:1:200-2000 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human skeletal muscle tissue showing brown staining of fibers.