

Product Datasheet

Mouse Cryopyrin ELISA Kit (orb408399)

Description	Mouse NACHT, LRR and PYD domains-containing protein 3(NLRP3) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	NLRP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8R4B8
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.