

## Product Datasheet

### Mouse CD274 ELISA Kit (orb408497)

**Description** Mouse Programmed cell death 1 ligand 1(CD274) ELISA kit

**Reactivity** Mouse

**Range** 31.25 pg/mL-2000 pg/mL

**Target** CD274

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [Q9EP73](#)

**Sensitivity** 7.81 pg/mL

**Expiration Date** Please enquire.