

## Product Datasheet

### Mouse PTGDS ELISA Kit (orb408505)

<b>Description</b>	Mouse Prostaglandin-H2 D-isomerase(PTGDS) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	4.7 ng/mL-300 ng/mL
<b>Target</b>	PTGDS
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">O09114</a>
<b>Sensitivity</b>	1.2 ng/mL
<b>Expiration Date</b>	Please enquire.