

Product Datasheet

Mouse PDIA4 ELISA Kit (orb408515)

Description	Mouse Protein disulfide-isomerase A4(PDIA4) ELISA kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	PDIA4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P08003
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.