

## Product Datasheet

### Mouse REG3A ELISA Kit (orb408573)

<b>Description</b>	Mouse Regenerating islet-derived protein 3-alpha(REG3A) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	3.12 ng/mL-200 ng/mL
<b>Target</b>	REG3A
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">O09037</a>
<b>Sensitivity</b>	0.78 ng/mL
<b>Expiration Date</b>	Please enquire.