

## Product Datasheet

### Mouse RSPO4 ELISA Kit (orb408592)

<b>Description</b>	Mouse R-spondin-4(RSPO4) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	34.4 pg/mL-2200 pg/mL
<b>Target</b>	RSPO4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q8BJ73</a>
<b>Sensitivity</b>	8.6 pg/mL
<b>Expiration Date</b>	Please enquire.